

Daily/ Weekly Test Sheet - Ultrasonic Cleaner (SHTM 01-05)

Tests to be carried out in accordance with Scottish Health Technical Memorandum (SHTM 01-05) Part B.

Table 1 - Record details

Hospital/ Location		Week beginning		Week number	
Department		Ref number		Serial number	

Daily Checks

Table 2 - Daily checks

Day	Filter and strainers removed and cleaned	Detergent levels	Drain chamber (Session 1)	Drain chamber (Session 2)
Monday				
Tuesday				
Wednesday				
Thursday				
Friday				
Saturday				
Sunday				

Weekly Checks

Table 3 - Weekly checks

Condition of door/ lid seal and interlocks:		Condition of bath, load holder/ basket:	
Protein residue test method used:		Batch number: Expiry date:	
Instrument(s) swabbed		Result	Pass/ Fail
Additional comments:			

Weekly Automatic Control Test (ACT) Check

Table 4 - Weekly ACT check. Note when not in use

Timer setting (min:sec)		Temperature setting (°C)					
Day	Volume of detergent (ml)	Water changed and degassed by	Cycle time (min:sec)	Indicated temperature (°C)	Satisfactory cleaning (visual check)	Fit for use name user	Fit for use signed user
Monday					Pass/ Fail		
Tuesday					Pass/ Fail		
Wednesday					Pass/ Fail		
Thursday					Pass/ Fail		
Friday					Pass/ Fail		
Saturday					Pass/ Fail		
Sunday					Pass/ Fail		

Detergent Details

Table 5 - Detergent details

Detergent	
Detergent name	
Detergent type	

Faults Record

Table 6 - Table to record faults (new or existing). Also enter in plant history record in FAC414-092

Fault comments