

## Daily/ Weekly Test Sheet - Small Steam Sterilizer with Type-N cycles (SHTM 01-05)

Tests to be carried out in accordance with Scottish Health Technical Memorandum (SHTM) 01-05 Part B.

Table 1 - Record details

Hospital/ Location		Week beginning		Week number	
Department		Ref number		Serial number	

### Automatic Control Test

Table 2 - Automatic control test - Note: SHTM 01-05 recommends a weekly test. Optional to update daily. SHTM 01-05 Part B recommends a production load of instruments is used for this test

Day	Cycle No	Sterilizing hold temp (°C)	Sterilizing hold pressure (bar/ kPa)	Sterilizing hold time	Result	Fit for use
Monday					Pass/ Fail	
Tuesday					Pass/ Fail	
Wednesday					Pass/ Fail	
Thursday					Pass/ Fail	
Friday					Pass/ Fail	
Saturday					Pass/ Fail	
Sunday					Pass/ Fail	

## Reservoir changes

Table 3 - Reservoir changes (where applicable). See SHTM 01-05 Part B 7.22 - using a small steam sterilizer which references Sterilized Water for Irrigation water quality used

Day	Cycle No when water changed	Comments	Water changed by
Monday			
Tuesday			
Wednesday			
Thursday			
Friday			
Saturday			
Sunday			

## Weekly safety checks

Table 4 - Table to record weekly safety checks

Door seal	Chamber/ shelves	Door closed interlock	Door pressure interlock
Satisfactory/ not satisfactory	Satisfactory/ not satisfactory	Satisfactory/ not satisfactory	Satisfactory/ not satisfactory

Faults record

Table 5 - Table to record faults (new or existing). Also enter in plant history record in FAC414-050.02

Fault comments

Test outcome

Table 6 - Table to record test outcome

Tested by	Date completed	Test result
		Satisfactory/ Not Satisfactory

